



Mouse Micro-Vent Rack Information

Standard Configurations: Opposite Side Plenums						
Size	Configuration*	Plenums	Width	Depth	Height**	Cage Size
30 Cage	6 W x 5 H SF	Opp	82.062" 2084	27.531" 699	58.500" 1486	10" x 19" x 6"
36 Cage	6 W x 6 H SF	Opp	82.062" 2084	27.531" 699	67.750" 1721	10" x 19" x 6"
70 Cage	7 W x 10 H SF	Rear	63.000" 1600	23.000" 584	80.062" 2034	7 1/2" x 11 " x 5 "
80 Cage	8 W x 10 H SF	Rear	71.500" 1816	23.000" 584	80.062" 2034	7 1/2" x 11 " x 5 "
84 Cage	7 W x 6 H DF	Opp	69.812" 1770	29.250" 743	52.312" 1329	7 1/2" x 11 " x 5 "
98 Cage	7 W x 7 H DF	Opp	69.812" 1770	29.250" 743	59.062" 1500	7 1/2" x 11 " x 5 "
112 Cage	7 W x 8 H DF	Opp	69.812" 1770	29.250" 743	65.812" 1672	7 1/2" x 11 " x 5 "
126 Cage	7 W x 9 H DF	Opp	69.812" 1770	29.250" 743	72.562" 1843	7 1/2" x 11 " x 5 "
140 Cage	7 W x 10 H DF	Opp	69.812" 1770	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "
160 Cage	8 W x 10 H DF	Opp	78.188" 1986	29.250" 743	79.812" 2027	7 1/2" x 11 " x 5 "

Optional Configurations: Same Side Plenums						
Size	Configuration*	Plenums	Width	Depth	Height**	Cage Size
84 Cage	7 W x 6 H DF	Same	71.031" 1804	29.250" 743	52.312" 1329	7 1/2" x 11 " x 5 "
98 Cage	7 W x 7 H DF	Same	71.031" 1804	29.250" 743	59.062" 1500	7 1/2" x 11 " x 5 "
112 Cage	7 W x 8 H DF	Same	71.031" 1804	29.250" 743	66.812" 1697	7 1/2" x 11 " x 5 "
126 Cage	7 W x 9 H DF	Same	71.031" 1804	29.250" 743	72.562" 1843	7 1/2" x 11 " x 5 "
140 Cage	7 W x 10 H DF	Same	71.031" 1804	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "
160 Cage	8 W x 10 H SF	Same	79.531" 2020	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "

*Rack Configuration: SF-Single Face, DF-Double Face

**Rack height: Add additional height of 10.50" (266.700mm) for rack mounted blower.

Mouse Micro-Vent Rack Information on following page.



Rat Micro-Vent Rack Information

Size	Configuration*	Plenums	Width	Depth	Height**	Cage Size
30 Cage	6 W x 5 H SF	Opp	82.000" 2083	27.500" 699	69.500" 1765	10 1/2" x 19" x 8"
36 Cage	6 W x 6 H SF	Opp	82.000" 2083	27.500" 699	80.000" 2032	10 1/2" x 19" x 8"
30 Cage	6 W x 5 H SF	Same	83.750" 2127	27.500" 699	69.000" 1753	10 1/2" x 19" x 8"
36 Cage	6 W x 6 H SF	Same	83.730" 2127	27.500" 699	80.250" 2038	10 1/2" x 19" x 8"
80 Cage	5 W x 8 H DF	Opp	72.812" 1849	36.000" 914	80.812" 2053	10" x 14" x 7"
80 Cage	5 W x 8 H DF	Same	75.562" 1920	36.000" 914	80.812" 2053	10" x 14" x 7"
40 Cage	5 W x 8 H SF	Rear	66.000" 1676	25.375" 645	80.812" 2053	10" x 14" x 7"

*Rack Configuration: SF-Single Face, DF-Double Face

**Rack height: Add additional height of 10.50" (266.700mm) for rack mounted blower.

Allentown reserves the right to change rack dimensions without notice.