

Micro-VENT System Caging All dimensions nominal	
<p> Model: MBS7115RHMV Mouse Cage: PC7115HT Wire Bar Lid: WBL7115SMD Reduced Height Micro-Barrier Top: MBT7115RH </p> <p> Overall dimension: 7 3/4"W x 12"D x 6 1/2"H Floor area: 67 sq.in. Shown with cage mounted valve </p>	
<p> Model: MBS75JRHMV Mouse Cage: "75 JAG" PC75JHT Wire Bar Lid: WBL7115SMD Reduced Height Micro-Barrier Top: MBT7115RH </p> <p> Overall dimension: 7 3/4"W x 12"D x 6 1/2"H Floor area: 75 square inches Shown with cage mounted valve </p>	
<p> Model: EBS2A75JFMV Mouse Cage: "75 JAG" PC75JHT Wire Bar Lid: WBLEBT External Bottle Assembly Micro-Barrier Top: EB27115F </p> <p> Overall dimension: 7 3/4"W x 12"D x 6 1/2"H Floor area: 75 square inches Shown with flap grommet for rack mounted valve (fixed or removable) </p>	
<p> Model: EBS2A75JFMV Mouse Cage: "75 JAG" PC75JHT Wire Bar Lid: WBLEBT External Bottle Assembly, Solid Top: EB27115ST </p> <p> Overall dimension: 7 3/4"W x 12"D x 6 1/2"H Floor area: 75 square inches Shown with flap grommet for rack mounted valve (fixed or removable) </p>	

Model: Mouse Large Cage: PC1019HT
Wire Bar Lid: WBL1019MSMD
Filter Micro-Barrier Top: MBT1019HT

Overall dimension: 10 3/4"W x 19 1/4"D x 8 7/8"H
Floor area: 153 square inches
Shown with cage mounted valve



Model: MBS10198MV
Rat Cage: PC10198HT
Wire Bar Lid: WBL1019RSMD
Filter Micro-Barrier Top: MBT1019HT

Overall dimension: 10 3/4"W x 19 1/4"D x 10 3/4"H
Floor area: 143 square inches
Shown with cage mounted valve



Model: MBS1014MV
Rat Cage: PC10147HT
Wire Bar Lid: WBL1014RSMD
Micro-Barrier Top: MBT1014HT

Overall dimension: 11.25" W x 15.5" D x 7.6"H
Floor area: 140 square inches
Shown with flap grommet for rack mounted valve
(fixed or removable)



P/NC System Caging

All dimensions nominal

Model: MBS75JSTMV
Mouse Cage: "75 JAG" PC75JHT
Wire bar lid: WBL7115SMD
Reduced Height Solid Micro-Barrier Top:
MBT7115RHSTPNCL

Overall dimension: 7 3/4"W x 12"D x 6 1/2"H
Floor area: 75 square inches
Shown with flap grommet for rack mounted valve
(fixed or removable)



Model: PCT2L12VCLMV
Mouse Cage Type II Long: PCT2LHT
External Bottle Assembly Micro-Barrier Top:
MBTTIICLPNC
Wire Bar Lid: WBLT2

Overall dimension: 8.4" W x 14.5" D x 7.3" H
Floor area: 80 square inches
Shown with flap grommet for rack mounted valve
(fixed or removable)




Model: PCT2L12VCLMV
Mouse Cage Type I Long: PCT2LHT
External Bottle Assembly Micro-Barrier Top:
MBTTIICLPNC
Wire Bar Lid: WBLT2

Overall dimension: 8.4" W x 14.5" D x 7.3" H
Floor area: 80 square inches
Shown with flap grommet for rack mounted valve
(fixed or removable)



Model: MBS1520MVST-PNC
Mouse Cage Multi-Species: PC1520HT
Hanging Feeders: F646
Solid Top: MBT1520ST-PNC

Overall dimension: 15"W x 19.5"D x 8.5"H
Floor areas: 210 square inches
Shown with flap grommet for rack mounted valve
(fixed or removable)



Allentown reserves the right to change rack dimensions without notice.