



## MiniDSR

Physical Characteristics	
Robot Type	Fanuc Arcmate 100 iBe, 6 kg.
Weight	138 kg. / 304.23 lb.
Installation Area	9 ft. <sup>2</sup>
Reach	1373 mm. / 54 in.
Power	
Required Power	480 v., 3-phase, 20 amp
Cage Capacity	
Quantity	1-2 cage per lift
Speed (1 Cage)	up to 300 cages per hour
Speed (2 Cage)	up to 600 cages per hour
Size	Wide variety
Related Equipment	
Disposal Station	Utilizes existing
Tunnel Washer	Utilizes existing
Cage Pallet	Required
Software	Included

Allentown reserves the right to change the dimensions without notice.