



RACK INFORMATION

inches/mm

Standard Configurations: Opposite Side Plenums						
Size	Configuration*	Plenums	Width	Depth	Height**	Cage Size
36 Cage	6 W x 6 H SF	Opp	82.000" 2084	27.500" 699	80.000" 2032	10 1/2" x 19" x 8"
70 Cage	7 W x 10 H SF	Rear	63.000" 1600	23.000" 584	80.062" 2034	7 1/2" x 11 " x 5 "
80 Cage	8 W x 10 H SF	Rear	71.500" 1986	23.000" 584	80.062" 2034	7 1/2" x 11 " x 5 "
84 Cage	7 W x 6 H DF	Opp	69.812" 1770	29.250" 743	52.312" 1329	7 1/2" x 11 " x 5 "
98 Cage	7 W x 7 H DF	Opp	69.812" 1770	29.250" 743	59.062" 1500	7 1/2" x 11 " x 5 "
112 Cage	7 W x 8 H DF	Opp	69.812" 1770	29.250" 743	65.812" 1672	7 1/2" x 11 " x 5 "
126 Cage	7 W x 9 H DF	Opp	69.812" 1770	29.250" 743	72.562" 1843	7 1/2" x 11 " x 5 "
140 Cage	7 W x 10 H DF	Opp	69.812" 1770	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "
160 Cage	8 W x 10 H DF	Opp	78.188" 1986	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "

Optional Configurations: Same Side Plenums						
Size	Configuration*	Plenums	Width	Depth	Height**	Cage Size
84 Cage	7 W x 6 H DF	Same	71.031" 1804	29.250" 743	52.312" 1329	7 1/2" x 11 " x 5 "
98 Cage	7 W x 7 H DF	Same	71.031" 1804	29.250" 743	59.062" 1500	7 1/2" x 11 " x 5 "
112 Cage	7 W x 8 H DF	Same	71.031" 1804	29.250" 743	66.812" 1697	7 1/2" x 11 " x 5 "
126 Cage	7 W x 9 H DF	Same	71.031" 1804	29.250" 743	72.562" 1843	7 1/2" x 11 " x 5 "
140 Cage	7 W x 10 H DF	Same	71.031" 1804	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "
160 Cage	8 W x 10 H DF	Same	79.531" 1804	29.250" 743	79.312" 2015	7 1/2" x 11 " x 5 "

*Rack Configuration: SF-Single Face, DF-Double Face

**Rack height: Add additional height of 10.50" (266.700mm) for rack mounted blower.

Allentown reserves the right to change rack dimensions without notice.