



TYPE I LONG RACK INFORMATION

Size	Configuration**	Plenums	Width	Depth	Height***
64 Cage	8 W x 8 H SF	Rear	1562 mm	659mm	1946mm
72 Cage	9 W x 8 H SF	Rear	1746mm	659mm	1946mm
80 Cage	10 W x 8 H SF	Rear	1930mm	659mm	1946mm
112 Cage	7 W x 8 H DF	Opp	1550mm	892mm	1946mm
128 Cage	8 W x 8 H DF	Opp	1734mm	892mm	1946mm
144 Cage	9 W x 8 H DF	Opp	1918mm	892mm	1946mm
160 Cage	10 W x 8 H DF	Opp	2102mm	892mm	1946mm

**Rack Configuration: SF-Single Face, DF-Double Face

***Rack height: Add additional height of 280mm for rack mounted blower.

Allentown reserves the right to change rack dimensions without notice.

TYPE II LONG RACK INFORMATION

Size	Configuration**	Plenums	Width	Depth	Height***
48 Cage	6 W x 8 H SF	Rear	59.00"/1499mm	25.688"/652mm	76.50"/1943mm
56 Cage	7 W x 8 H SF	Rear	68.25"/1734mm	25.688"/652mm	76.50"/1943mm
64 Cage	8 W x 8 H SF	Rear	77.50"/1969mm	25.688"/652mm	76.50"/1943mm
96 Cage	6 W x 8 H DF	Opp	65.81"/1672mm	34.625"/879mm	76.50"/1943mm
112 Cage	7 W x 8 H DF	Opp	75.06"/1907mm	34.625"/879mm	76.50"/1943mm
128 Cage	8 W x 8 H DF	Opp	84.31"/2142mm	34.625"/879mm	76.50"/1943mm
144 Cage*	8 W x 9 H DF	Opp	84.31"/2142mm	34.625"/879mm	80.94"/2056mm
96 Cage	6 W x 8 H DF	Same	66.66"/1693mm	34.625"/879mm	76.50"/1943mm
112 Cage	7 W x 8 H DF	Same	75.90"/1928mm	34.625"/879mm	76.50"/1943mm
128 Cage	8 W x 8 H DF	Same	85.16"/2163mm	34.625"/879mm	76.50"/1943mm
144 Cage*	8 W x 9 H DF	Same	85.16"/2163mm	34.625"/879mm	80.94"/2056mm

*Only available with 127mm high cage

**Rack Configuration: SF-Single Face, DF-Double Face

***Rack height: Add additional height of 11"/280mm for rack mounted blower.

Allentown reserves the right to change rack dimensions without notice.

Type III and Type III High Information on following page.



TYPE III & TYPE III HIGH RACK INFORMATION

Size	Configuration	Length	Width	Height
24 Cage	4 W x 6 H	691mm	1295mm	1938mm
30 Cage	5 W x 6 H	691mm	1594mm	1938mm
36 Cage	6 W x 6 H	691mm	1892mm	1938mm

Allentown reserves the right to change rack dimensions without notice.