



SINGLE FACE WITH REAR PLENUM

Rack Configurations	SF/DF	Number of Cages	Length IN/mm	Width IN/mm	Height IN/mm
6W x 7H	SF	42	52.250 / 1327	26.000 / 660	60.897 / 1547
6W X 8H	SF	48	52.250 / 1327	26.000 / 660	68.272 / 1734
6W X 9H	SF	54	52.250 / 1327	26.000 / 660	75.647 / 1921
6W X 10H	SF	60	52.250 / 1327	26.000 / 660	83.022 / 2109
7W x 7H	SF	49	60.500 / 1537	26.000 / 660	60.897 / 1547
7W X 8H	SF	56	60.500 / 1537	26.000 / 660	68.272 / 1734
7W X 9H	SF	63	60.500 / 1537	26.000 / 660	75.647 / 1921
7W X 10H	SF	70	60.500 / 1537	26.000 / 660	83.022 / 2109
8W x 7H	SF	56	68.750 / 1746	26.000 / 660	60.897 / 1547
8W X 8H	SF	64	68.750 / 1746	26.000 / 660	68.272 / 1734
8W X 9H	SF	72	68.750 / 1746	26.000 / 660	75.647 / 1921
8W X 10H	SF	80	68.750 / 1746	26.000 / 660	83.022 / 2109
9W x 7H	SF	63	77.000 / 1956	26.000 / 660	60.897 / 1547
9W X 8H	SF	72	77.000 / 1956	26.000 / 660	68.272 / 1734
9W X 9H	SF	81	77.000 / 1956	26.000 / 660	75.647 / 1921
9W X 10H	SF	90	77.000 / 1956	26.000 / 660	83.022 / 2109

DOUBLE FACE WITH SAME SIDE PLENUM

Rack Configurations	SF/DF	Number of Cages	Length IN/mm	Width IN/mm	Height IN/mm
6W x 7H	DF	84	56.172 / 1427	35.500 / 902	60.897 / 1547
6W X 8H	DF	96	56.172 / 1427	35.500 / 902	68.272 / 1734
6W X 9H	DF	108	56.172 / 1427	35.500 / 902	75.647 / 1921
6W X 10H	DF	120	56.172 / 1427	35.500 / 902	83.022 / 2109
7W x 7H	DF	98	64.422 / 1636	35.500 / 902	60.897 / 1547
7W X 8H	DF	112	64.422 / 1636	35.500 / 902	68.272 / 1734
7W X 9H	DF	126	64.422 / 1636	35.500 / 902	75.647 / 1921
7W X 10H	DF	140	64.422 / 1636	35.500 / 902	83.022 / 2109

8W x 7H	DF	112	72.672 / 1846	35.500 / 902	60.897 / 1547
8W X 8H	DF	128	72.672 / 1846	35.500 / 902	68.272 / 1734
8W X 9H	DF	144	72.672 / 1846	35.500 / 902	75.647 / 1921
8W X 10H	DF	160	72.672 / 1846	35.500 / 902	83.022 / 2109
9W x 7H	DF	126	80.922 / 2055	35.500 / 902	60.897 / 1547
9W X 8H	DF	144	80.922 / 2055	35.500 / 902	68.272 / 1734
9W X 9H	DF	162	80.922 / 2055	35.500 / 902	75.647 / 1921
9W X 10H	DF	180	80.922 / 2055	35.500 / 902	83.022 / 2109

Allentown reserves the right to change rack dimensions without notice.